




Presents the new 3-Dimensional
Open Loop Designed

AndraSnare

&

AndraSnare Micro


AndraSnare with Open Loop Design



- Real 3-Dimensional retrievalability; to provide best manipulation options!
- High upstanding hoop force for excellent procedural success !


* Closed Loop Designed Snares can not manipulate around a foreign body!

AndraSnare with preshaped Guiding Catheter




- Multi Purpose pre-shaped Guiding catheter
Increases 3-Dimensional captureability and eases manipulation into target vessel. For high procedural success
- The BIGG Torque device allows the physician to easily manipulate in all 3-Dimensional directions

AndraSnare with advanced Accessories



- SLIGHT insertion tool to easy load the snare and minimize blood loss during intervention increasing safety of the procedure
- Long supportive sleeve to provide high stability during manipulation & minimizing cost! (Especially when a new catheter is needed due to breakage.)

AndraSnare with High Product Variability and extra long UBD




- AndraSnare comes in the highest product variability!
- AndraSnare has best Sheath compatibilities! For better patient care.
- AndraSnare has longest UBD (4 years) to minimize cost!

AndraSnare & Retrieval Techniques

- Choose loop size of the Andrasnare smaller than vessel diameter to allow full expansion

Proximal Grab Technique



stent

Distal Wire Grab Technique

AndraSnare Other Retrieval Techniques

advance catheter in standard technique close to target



exchange guide wire to AndraSnare



advance / rotate the AndraSnare to catch the foreign body



Technical Description

Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
AS-5	5mm (3-5mm vessel)	125	110	4
AS-10	10mm (6-10mm vessel)	125	110	4
AS-15	15mm (11-15mm vessel)	125	110	4
AS-20	20mm (16-20mm vessel)	125	110	4
AS-25	25mm (21-25mm vessel)	125	110	5
AS-30	30mm (26-30mm vessel)	125	110	5
AS-35	35mm (31-35mm vessel)	125	110	5

Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
AS-25/65	25mm (21-25mm vessel)	65	50	5
AS-25/150	25mm (21-25mm vessel)	150	135	5

Product Code	Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
ASM-2	2mm (1-2mm vessel)	175	150	3-2.3
ASM-4	4mm (3-4mm vessel)	175	150	3-2.3
ASM-7	7mm (5-7mm vessel)	175	150	3-2.3

