




Presents the new CoCr (Cobalt  
Chromium) Hybrid Cell Design  
XL & XXL Stents

# AndraStent XL & XXL with Hybrid Cell Design



- Real Hybrid Cell design combines
  - flexibility of Open Cell Design and
  - radial strength of Closed Cell Design Stents!
- Optimal Scaffolding optimizes the clinical result!

\*Hybrid Cell Design allows full flexibility pre,- and post deployment!

# AndraStent XL & XXL Advantages of Cobalt Chromium (CoCr) Stent Material!




- AndraStent XL & XXL (CoCr) reduces sheath compatibility

- AndraStent XL & XXL (CoCr) maximizes radial strength

- Cobalt reduces implant material



# AndraStent XL & XXL broad Product variability!



- AndraStent XL  
Expansion Range  
< 25 mm
- AndraStent XXL  
Expansion Range  
< 32 mm
- Available Lengths  
- 13, 17, 21, 26, 30  
39, 48, 57 mm

# Technical Description

- Cobalt Chromium (CoCr) MP35N Stent Material
- Lowest Stent Profile
- Combined With High Radial Strength And Flexibility

Order Codes AndraStent XL and XXL									
	Stent Length								
	13 mm	17 mm	21 mm	26 mm	30 mm	39 mm	43 mm	48 mm	57 mm
Andrastent XL	AS13XL	AS17XL	AS21XL	AS26XL	AS30XL	AS39XL	AS43XL	AS48XL	AS57XL
Andrastent XXL	-----	AS17XXL	AS21XXL	AS26XXL	AS30XXL	AS39XXL	AS43XXL	AS48XXL	AS57XXL