


**AndraValvulotome<sup>®</sup>**  
& Andravalvulotome<sup>®</sup> OTW

Nitinol System (4 and 5 F)

# AndraValvulotome®




- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves
- Atraumatic tip

## Order Codes AndraValvulotome®

Catalogue Number	Max. Diameter Basket [mm]	Max. Diameter Blades [mm]	Diameter Catheter [mm]	Usable Length [mm]	Guide Wire
AVT-Tip-5-10/1000	8	4.2	1.67	1000	-

Packaging Unit: 1

# AndraValvulotome® OTW



- The OTW-system permits precise guidance of the valvulotome
- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves

## Order Codes AndraValvulotome® OTW

Catalogue Number	Max. Diameter Basket [mm]	Max. Diameter Blades [mm]	Diameter Catheter [mm]	Usable Length [mm]	Guide Wire
AVT-5F-8/1000	8	6.8	1.67	1000	0.018 "
AVT-4F-8/900	6	5.5	1.37	900	0.018 "

Packaging Unit: 1

## **Andramed GmbH**

Schiesswiesenstr. 18  
72766 Reutlingen, Germany  
Tel.: +49 (0) 7127 811 57 15  
Fax: +49 (0) 7127 811 57 11  
[info@andramed.com](mailto:info@andramed.com)  
[www.andramed.com](http://www.andramed.com)



Made in Germany  
Rev.07/2023-10-25

