


AndraValvulotome[®]
& Andravalvulotome[®] OTW

Nitinol system (4 and 5 Fr)

AndraValvulotome®




- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves
- Atraumatic tip

Order Codes AndraValvulotome®

Order Code	Max. Hoop Ø in mm	Blade Range Ø in mm	Protective Sheath Outer Ø	Shaft Length in mm
AVT-TIP-5F-10/1000	6	3 – 6	5 Fr/1.67 mm	1000

Packaging Unit: 1

AndraValvulotome® OTW



- The OTW-system permits precise guidance of the valvulotome
- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves

Order Codes AndraValvulotome® OTW

Order Code	Max. Hoop Ø in mm	Blade Range Ø in mm	Protective Sheath Outer Ø	Shaft Length in mm	Guide Wire
AVT-5F-8/1000	8	3 – 8	5 Fr/1.67	1000	0.018 "
AVT-4F-8/900	4.5	2 – 4	4 Fr/1.33	900	0.014 "

Packaging Unit: 1

Andramed GmbH

Schiesswiesenstr. 18
72766 Reutlingen, Germany
Tel.: +49 (0) 7127 811 57 15
Fax: +49 (0) 7127 811 57 11
info@andramed.com
www.andramed.com

Made in Germany
Rev. 05/2020

