


**AndraValvulotome<sup>®</sup>**  
& AndraValvulotome<sup>®</sup> OTW

Nitinol System (4 and 5 Fr)

# AndraValvulotome®




- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves
- Atraumatic tip

## Order Codes AndraValvulotome®

| Catalogue Number   | Max. Diameter Basket [mm] | Max. Diameter Blades [mm] | Diameter Catheter [mm] | Usable Length [mm] | Guide Wire |
|--------------------|---------------------------|---------------------------|------------------------|--------------------|------------|
| AVT-Tip-5F-10/1000 | 8                         | 4.2                       | 1.67                   | 1000               | -          |

Packaging Unit: 1

# AndraValvulotome® OTW



- The OTW-system permits precise guidance of the valvulotome
- Automatic self-sizing and self-centering technology because of innovative basket geometry and material
- The automatic adjustment to the vein diameter ensures cutting of the various valves

## Order Codes AndraValvulotome® OTW

| Catalogue Number | Max. Diameter Basket [mm] | Max. Diameter Blades [mm] | Diameter Catheter [mm] | Usable Length [mm] | Guide Wire |
|------------------|---------------------------|---------------------------|------------------------|--------------------|------------|
| AVT-5F-8/1000    | 8                         | 6.8                       | 1.67                   | 1000               | 0.018 "    |
| AVT-4F-8/900     | 6                         | 5.5                       | 1.37                   | 900                | 0.018 "    |

Packaging Unit: 1

## **Andramed GmbH**

Schiesswiesenstr. 18  
72766 Reutlingen, Germany  
Tel.: +49 (0) 7127 811 57 15  
Fax: +49 (0) 7127 811 57 11  
[info@andramed.com](mailto:info@andramed.com)  
[www.andramed.com](http://www.andramed.com)



Made in Germany  
Rev.06/2023-03-06

